

Viscosity & Density Testing

MODELS

KVD-70Xe



Purpose-Built for Viscosity & Density

KVD-70Xe is a dedicated analyzer engineered to provide highly precise viscosity and density measurement for jet fuels and diesel, supporting laboratories that demand accuracy, stability, and confidence in their results.

- **Method-Driven Performance**
- **Designed for Jet fuel & Diesel Laboratories**
- **Compact Footprint for Efficient Integration**

Included Functions	KVD-70Xe
D7945 Viscosity @ -20°C & -40°C	✓
D7945 Viscosity @ 40°C	✓
D7777 Density	✓
Renewable Diesel (HVO)	✓
Sustainable Aviation Fuel (SAF)	✓
Sample Introduction	Automatic Single Injection
Optional	30-Position Autosampler